



**miR-Direct®: assay for miRNA
quantification from human plasma
and blood**

Experimental Protocol

support@somagenics.com
www.somagenics.com

Table of Contents

I.	Overview	5
II.	miR-Direct® Kit Contents	6
III.	User Supplied Reagents and Laboratory Equipment (not included)	7
IV.	Experimental Protocol	8
	Subsection A: Plasma/Serum Protocol	9
	Subsection B: Whole Blood Protocol	14
V.	Appendix	19

For Research Use Only.

Information in this document is subject to change without notice.

SOMAGENICS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL SOMAGENICS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT SOMAGENICS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

NOTICES TO PURCHASER:

LIMITED LICENSE

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of SomaGenics, Inc. For information on obtaining additional rights, please contact: bd@somagenics.com.

miR-ID®/miR-Direct® CUSTOMER USE AGREEMENT

When purchasing products incorporating SOMAGENICS' miR-ID® and miR-Direct® technologies you agree to the following terms and conditions. These terms and conditions constitute the use agreement ("Use Agreement") for products and materials containing miR-ID® and miR-Direct® sample preparation products, primer designs, primer sequences, primers, probe designs, probes and/or assays (together "miR-ID® and miR-Direct® Materials"). This Use Agreement creates a binding contract between you and SOMAGENICS for the purchase and use of our miR-ID® and miR-Direct® Materials.

Under this Use Agreement, SOMAGENICS retains all right, title and interest in and to the miR-ID® and miR-Direct® Materials. All custom services derived from miR-ID® and miR-Direct® Materials are provided on a non-exclusive basis (see www.somagenics.com; Terms and Conditions for Custom Services). miR-ID® and miR-Direct® Materials are provided for Research Use Only and no other rights are conveyed by this Use Agreement. You are expressly prohibited from independently recreating miR-ID® and miR-Direct® Materials which are proprietary to SOMAGENICS. If you require additional rights please contact: bd@somagenics.com.

This product is covered by issued and pending patents.

TRADEMARKS

miR-ID® and miR-Direct® are the Registered Trademarks of SomaGenics, Inc. All other brands and names contained herein are the property of their respective owners.


support@somagenics.com

www.somagenics.com

miR-Direct[®]: miRNA Detection from Human Plasma and Blood

I. Overview

The miR-Direct[®] assay is designed to detect miRNAs present in plasma and blood (50 to 400 μ L). The first part of the assay is performed in one tube and entails miRNA release from plasma or blood, miRNA capture, washing steps, and miRNA release with simultaneous circularization. The second part is miRNA quantification by SomaGenics' miR-ID[®] technology, consisting of reverse transcription by rolling circle amplification and real-time qPCR using SYBR green. All steps up to quantitative qPCR can be performed in multiplex for all miRNAs of interest.



support@somagenics.com
www.somagenics.com

II. miR-Direct[®] Kit Contents

Store at ambient temperature

1. Buffer 1*

Store at 4°C

1. Wash buffer
2. Buffer 3

Store at -20°C

1. miR-cel-39 synthetic spike-in miRNA control
2. Spike-in Dilution Mix
3. Protease Solution**
4. Buffer 2
5. Probe Mix
6. Circularization Reaction Mix
7. RT Reaction Mix
8. RT Primers (or customized RT primer mix for all miRNAs of interest)
9. qPCR Primer Pairs (miRNA-specific)
10. Nuclease-free Water

* CAUTION: Irritant. Contains guanidine hydrochloride. Use proper personal protective equipment. Do not use bleach.

** Keep at -20°C upon receipt. Store at 4°C after thawing for first use.

III. User-supplied Reagents, Consumables, and Laboratory Equipment (not included)

Reagents

- Streptavidin Magnetic beads (New England BioLabs, cat# S1420S)
- SuperScript II Reverse Transcriptase (ThermoFisher Scientific, cat#18064022)
- CircLigase II (Epicentre (Illumina) cat# CL9021K)
- RNaseOUT (ThermoFisher Scientific, cat# 10777-019)
- 2X iQ SYBR Green Supermix (BioRad, cat#170-8880)

Consumables

- Pipettes and nuclease-free filter pipette tips for volumes 2 to 100 μ L.
- RNase-free microfuge tubes (1.5 μ L).
- PCR tubes or PCR strip tubes.
- Micro-Amp Fast Optical 96-Well Reaction Plate, for example ThermoFisherScientific # 4346907
- Micro-Amp Clear Adhesive Film, for example ThermoFisherScientific # 4311971

Equipment

- Bench top tube rack for 1.5 mL tubes
- Bench top PCR plate rack
- Ice bucket (and ice)
- 12-tube Magnetic Separation Rack, for example NEB cat# S1509S
- Centrifuge for 1.5 mL and PCR tubes
- Centrifuge for 96-well PCR plates
- End-over-end rotator (at 4°C)
- Water bath/heat block (25°C)
- Water bath/heat block (37°C)
- Water bath/heat block (60°C)
- Thermal Cycler
- RT-PCR Thermal Cycler, for example ABI 7500 Fast System

IV. Experimental Protocol

1. Component preparation

Resuspend magnetic beads by vortexing them at high speed and then placing them on a tube rotator at 4°C for at least 120 minutes.

NOTE: Magnetic beads tend to aggregate while in storage, so it is imperative that they are fully dispersed in solution before they are used.

1.1. Place at room temperature to thaw:

- miR-cel-39 spike-in miRNA
- Spike-in Dilution Mix
- Protease Solution* (first use only -- **do not vortex**)
- Buffer 2
- Probe Mix
- Circularization Reaction Mix
- RT Reaction Mix
- RT Primer Mix
- Nuclease-free Water
- Wash Buffer (although not frozen, wash buffer should be at room temperature prior to use)

Protease solution is shipped frozen. Store at -20°C until first use. After thawing for the first time, store at 4°C (DO NOT REFREEZE**).*

After thawing, all reagents (excluding wash buffer) should be kept on ice while not in use. Store enzymes at -20°C until they are needed, except protease solution, which should be kept at 4°C

Return all reagents and enzymes to their proper storage conditions after use.

1.2. Store all enzymes at -20°C (except protease solution) until they are needed in the protocol.

IV Subsection A: Plasma / Serum Protocol

2. Thaw plasma samples and prepare Spike-in Mix

2.1. Thaw plasma samples at 25°C. Use plasma IMMEDIATELY after it is thawed.

2.2. Preparation of Spike-in Mix

While plasma is thawing, prepare the Spike-in Mix according to the table below for the volume of plasma input into the assay. To prepare a master mix to analyze multiple samples, scale the volumes listed below to the desired number of plasma samples to be analyzed. Vortex the Spike-in Mix for 3 seconds after the addition of the components.

	Plasma Volume		
	50 µL	200 µL	400 µL
miR-cel-39	4 µL	4 µL	4 µL
Spike-in Dilution Mix	1.2 µL	4.8 µL	9.6 µL

3. Plasma lysis

NOTE: Work *one tube at a time*.

3.1. Depending on the amount of plasma used, add the proper reagent volume to the sample in the order listed below. Vortex samples for 3 seconds after the addition of each reagent.

	Plasma Volume		
	50 µL	200 µL	400 µL
Buffer 1*	50 µL	200 µL	400 µL
Protease Solution	2 µL	8 µL	16 µL
Spike-in Mix	5.2 µL	8.8 µL	13.6 µL

* CAUTION: Irritant. Contains guanidine hydrochloride.

3.2. Incubate samples at 25°C for 60 min

4. miRNA capture

NOTE: Work *one tube at a time*.

4.1. Depending on the amount of plasma used, add the proper reagent volume to the sample in the order listed below. Vortex samples for 3 seconds after the addition of each reagent.

support@somagenics.com

www.somagenics.com

	Plasma Volume		
	50 μL	200 μL	400 μL
Buffer 2	11.5 μL	46 μL	92 μL
Probe Mix	6.4 μL	6.4 μL	6.4 μL

- 4.2. Incubate at 37°C for 60 min
- 4.3. Add 20 μL magnetic beads (and Buffer 3 if indicated) and vortex each tube for 3 seconds after addition.

	Plasma Volume		
	50 μL	200 μL	400 μL
Magnetic Beads	20 μL	20 μL	20 μL
Buffer 3	275 μL	none	none

- 4.4. Incubate at room temperature for 15 minutes
- Pulse spin tubes at 325 x g (approximately 3 seconds) to collect the liquid. **Do not allow beads to pellet.**

5. Washing steps

- 5.1. Place all samples on the magnetic rack for 10 minutes.

NOTE: Work *one tube at a time*.

- 5.2. While keeping the tube on the rack, carefully remove the solution from the beads by pipetting from the top of the tube, keeping the pipette tip against the tube wall opposite of the magnetic bead pellet. Be very careful not to disturb the beads.

Discard supernatant.

- 5.3. Remove the tube from the magnetic rack.
- 5.4. Add 500 μL room temperature wash buffer and vortex for 3 seconds. Leave tube in regular tube rack.

Repeat steps 5.1 to 5.4 two more times for each tube except that only 470 μL are removed after the final wash (leaving 30 μL buffer with the magnetic beads).

Keep tubes on a bench top rack until all samples are processed.

- 5.5. Pulse vortex tubes to resuspend the magnetic bead pellet.

support@somagenics.com
www.somagenics.com

- 5.6. Pulse spin tubes at 325 x g (approximately 3 seconds) to collect the liquid. **Do not allow beads to pellet.**

Leave samples on benchtop rack.

6. miRNA circularization

- 6.1. Mix the following reagents in a separate tube (same volume for all plasma input volumes). When assaying multiple samples, make a combined master mix and allow 10% overage for pipetting loss. Gently mix components and pulse-spin to collect liquid on bottom. Keep on ice.

Reagent	Volume
Circularization Reaction Mix	5 μ L
RNaseOUT	0.6 μ L
CircLigase II Enzyme	0.4 μ L
Nuclease-free Water	4 μ L
Total	10.0 μL

NOTE: Work *one tube at a time*.

- 6.2. Place sample tube from step 5.6. on the magnetic rack. When beads have aggregated, solution will be clear.
- 6.3. While keeping the tube on the rack, carefully remove the solution from the beads, keeping the pipette tip against the tube wall opposite of the magnetic bead pellet. Be very careful not to disturb the beads.
- 6.4. Immediately after removing wash buffer, remove the tube from the magnetic rack and pipet 10 μ L circularization mixture from step 6.1. onto the bead pellet. Pipette up and down 10 times to fully resuspend the beads into the circularization mixture, making sure that there are no beads left on the side of the tube.
- 6.5. Repeat steps 6.2 to 6.4 with the remaining samples.
- 6.6. Pulse spin tubes at 325 x g (approximately 3 seconds) to collect the liquid. **Do not allow beads to pellet.**
- 6.7. Incubate beads in circularization mixture at 60°C for 15 minutes.

7. Reverse transcription

- 7.1. Prepare the RT master mix while the circularization reaction is proceeding: mix the following reagents in a separate tube. When assaying multiple samples, make a combined master mix and allow 10% overage for pipetting loss. Gently mix components and pulse-spin to collect liquid on bottom. Keep on ice.

Reagent	Volume
RT Reaction Mix	25.6 μ L
RT Primer Mix	32 μ L
SuperScript II Enzyme	1.6 μ L
Nuclease-free Water	10.8 μ L
Total	70.0 μL

- 7.2. Remove tubes from 60° (from step 6.7.) and place them on the magnetic rack.

NOTE: Work *one tube at a time*.

- 7.3. After magnetic beads have migrated to the tube wall, carefully pipet the supernatant (10 μ L; free of magnetic beads) into new, labeled PCR tubes.
- 7.4. Add 70 μ L of RT master mix to each sample (80 μ L final volume).
- 7.5. Mix each sample gently and collect liquid by briefly spinning the tubes at 325 x g.
- 7.6. Place sample PCR tubes into thermal cycler and run the following parameter values:

Step type	Temperature	Time
HOLD	42°C	60 min
HOLD	75°C	15 min
HOLD	4°C	∞

***Stopping Point*:** RT reactions may be used directly in qPCR or stored at -20°C.

8. Quantitative RT PCR

- 8.1. If RT products were frozen, thaw samples and pulse vortex to mix and pulse spin to collect liquid.
- 8.2. Thaw 2x iQ SYBR green Master Mix and ROX dyes before use.

support@somagenics.com
www.somagenics.com

NOTE: SYBR green and ROX dyes are **light sensitive**. Perform all steps involving SYBR Green and ROX dyes away from direct light sources and keep tubes in dark when not in use.

8.3. Mix qPCR reagents as listed below.

A separate master mix needs to be prepared for each assayed miRNA. Calculate master mix components for each miRNA by multiplying the reagent volume by the number of samples and by the number of technical repeats (triplicates are recommended for quantitative RT-PCR assays). Allow 10% overage for pipetting loss. Mix gently and spin down the qPCR master mix.

Reagent	Volume
2X iQ SYBR Green Supermix	10 μ L
PCR Primers (miRNA-specific)	0.6 μ L
qPCR Reference Dye (ROX)	0.3 μ L
Nuclease-free Water	7.1 μ L
Total volume	18 μL

8.4. Pipet 18 μ L of the appropriate master mix into the well of a qPCR plate.

8.5. Pipet 2 μ L cDNA from the RT reaction (after 7.6. is completed) into the appropriate wells of the qPCR plate.

8.6. Seal qPCR plate with optical seal.

8.7. Centrifuge plate at 4000 rpm for 20-30 seconds.

8.8. Place plate in the ABI 7500 thermal cycler. Run qPCR using the following profile:

Step Type	Temperature	Time
HOLD	50°C	2 min
HOLD	95°C	3 min
CYCLE (40 cycles)	95°C	15 sec
	56°C	45 sec
	66°C + 0.1 per cycle	30 sec
	67°C + 0.2 per cycle*	30 sec
Include standard melting curve profile at the end of 40 cycles		

IV Subsection B: Whole Blood Protocol

2. Prepare Spike-in Mix

2.1. Preparation of Spike-in Mix

Prepare the Spike-in Mix according to the table below. To prepare a master mix to analyze multiple samples, scale the volumes listed below to the desired number of plasma samples to be analyzed. Vortex the Spike-in Mix for 3 seconds after the addition of the components.

Component	Volume to add
miR-cel-39	4 μL
Spike-in Dilution Mix	9.6 μL

3. Blood lysis

NOTE: Work *one tube at a time*.

3.1. Depending on the amount of blood used, add the proper reagent volume to the sample in the order listed below. Vortex samples for 3 seconds after the addition of each reagent.

Component	Plasma Volume	
	25 μL	200 μL
Buffer 1*	400 μL	200 μL
Protease Solution	16 μL	8 μL
Spike-in Mix	13.6 μL	8.8 μL
Buffer 3	345 μL	370 μL

* CAUTION: Irritant. Contains guanidine hydrochloride.

3.2. Incubate samples at 25°C for 60 min

4. miRNA capture

NOTE: Work *one tube at a time*.

4.1. Add the proper reagent volume to the sample in the order listed below. Vortex samples for 3 seconds after the addition of each reagent.

Component	Volume
Buffer 2	92 μL
Spike-in Dilution Mix	6.4 μL

- 4.2. Incubate at 37°C for 60 min
- 4.3. Add 20 µL magnetic beads and vortex each tube for 3 seconds after addition.
- 4.4. Incubate at room temperature for 15 minutes

Pulse spin tubes at 325 x g (approximately 3 seconds) to collect the liquid. **Do not allow beads to pellet.**

5. Washing steps

- 5.1. Place all samples on the magnetic rack for 10 minutes.

NOTE: Work *one tube at a time*.

- 5.2. While keeping the tube on the rack, carefully remove the solution from the beads, keeping the pipette tip against the tube wall opposite of the magnetic bead pellet. Be very careful not to disturb the beads.

Discard supernatant.

- 5.3. Remove the tube from the magnetic rack.
- 5.4. Add 500 µL room temperature wash buffer and vortex for 3 seconds. Leave tube in regular tube rack.

Repeat steps 5.1 to 5.4 two more times for each tube except that only 470 µL are removed after the final wash (leaving 30 µL buffer with the magnetic beads).

Keep tubes on a bench top rack until all samples are processed.

- 5.5. Pulse vortex tubes to resuspend the magnetic bead pellet.
- 5.6. Pulse spin tubes at 325 x g (approximately 3 seconds) to collect the liquid. **Do not allow beads to pellet.**

Leave samples on benchtop rack.

6. miRNA circularization

- 6.1. Mix the following reagents in a separate tube (same volume for all blood input volumes). When assaying multiple samples, make a combined master mix and allow 10% overage for pipetting loss. Gently mix components and pulse-spin to collect liquid on bottom. Keep on ice.

support@somagenics.com
www.somagenics.com

Reagent	Volume
Circularization Reaction Mix	5 μ L
RNAseOUT	0.6 μ L
CircLigase II Enzyme	0.4 μ L
Nuclease-free Water	4 μ L
Total	10.0 μL

NOTE: Work *one tube at a time*.

- 6.2. Place sample tube from step 5.6. on the magnetic rack. When beads have aggregated, solution will be clear.
- 6.3. While keeping the tube on the rack, carefully remove the solution from the beads, keeping the pipette tip against the tube wall opposite of the magnetic bead pellet. Be very careful not to disturb the beads.
- 6.4. Immediately after removing wash buffer, remove the tube from the magnetic rack and pipet 10 μ L circularization mixture from step 6.1. onto the bead pellet. Pipette up and down 10 times to fully resuspend the beads into the circularization mixture, making sure that there are no beads left on the side of the tube.
- 6.5. Repeat steps 6.2 to 6.4 with the remaining samples.
- 6.6. Pulse spin tubes at 325 x g (approximately 3 seconds) to collect the liquid. **Do not allow beads to pellet.**
- 6.7. Incubate beads in circularization mixture at 60°C for 15 minutes.

7. Reverse transcription

- 7.1. Prepare the RT master mix while the circularization reaction is proceeding: mix the following reagents in a separate tube. When assaying multiple samples, make a combined master mix and allow 10% overage for pipetting loss. Gently mix components and pulse-spin to collect liquid on bottom. Keep on ice.

Reagent	Volume
RT Reaction Mix	25.6 μ L
RT Primer Mix	32 μ L
SuperScript II Enzyme	1.6 μ L
Nuclease-free Water	10.8 μ L
Total	70.0 μL

7.2. Remove tubes from 60° (from step 6.7.) and place them on the magnetic rack.

NOTE: Work *one tube at a time*.

7.3. After magnetic beads have migrated to the tube wall, carefully pipet the supernatant (10 μ L; free of magnetic beads) into new, labeled PCR tubes.

7.4. Add 70 μ L of RT master mix to each sample (80 μ L final volume).

7.5. Mix each sample gently and collect liquid by briefly spinning the tubes at 325 x g.

7.6. Place sample PCR tubes into thermal cycler and run the following parameter values:

Step type	Temperature	Time
HOLD	42°C	60 min
HOLD	75°C	15 min
HOLD	4°C	∞

NOTE (optional atopping point): RT reactions may be used directly in qPCR or stored at -20°C.

8. Quantitative RT PCR

8.1. If RT products were frozen, thaw samples and pulse vortex to mix and pulse spin to collect liquid.

8.2. Thaw 2x iQ SYBR green Master Mix and ROX dyes before use.

NOTE: SYBR green and ROX dyes are **light sensitive**. Perform all steps involving SYBR Green and ROX dyes away from direct light sources and keep tubes in dark when not in use.

8.3. Mix qPCR reagents as listed on the following page.


A separate master mix needs to be prepared for each assayed miRNA. Calculate master mix components for each miRNA by multiplying the reagent volume by the number of samples and by the number of technical repeats (triplicates are recommended for quantitative RT-PCR assays). Allow 10% overage for pipetting loss. Mix gently and spin down the qPCR master mix.

Reagent	Volume
2X iQ SYBR Green Supermix	10 μ L
PCR Primers (miRNA-specific)	0.6 μ L
qPCR Reference Dye (ROX)	0.3 μ L
Nuclease-free Water	7.1 μ L
Total volume	18 μL

- 8.4. Pipet 18 μ L of the appropriate master mix into the well of a qPCR plate.
- 8.5. Pipet 2 μ L cDNA from the RT reaction (after 7.6. is completed) into the appropriate wells of the qPCR plate.
- 8.6. Seal qPCR plate with optical seal.
- 8.7. Centrifuge plate at 4000 rpm for 20-30 seconds.
- 8.8. Place plate in the ABI 7500 thermal cycler. Run qPCR using the following profile:

Step Type	Temperature	Time
HOLD	50°C	2 min
HOLD	95°C	3 min
CYCLE (40 cycles)	95°C	15 sec
	56°C	45 sec
	66°C + 0.1 per cycle	30 sec
	67°C + 0.2 per cycle*	30 sec
Include standard melting curve profile at the end of 40 cycles		

V. Appendix: Overview of Working Steps



support@somagenics.com
www.somagenics.com